

PACO47174

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

May act as a regulator of proteasome.

**Gene ID:**

CYREN

**Uniprot**

Q9BWK5

**Synonyms:**

Cell cycle regulator of non-homologous end joining (Cell cycle regulator of NHEJ) (Modulator of retrovirus infection homolog), CYREN, C7orf49 MRI

**Immunogen:**


Recombinant Human Cell cycle regulator of non-homologous end joining protein (1-157AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human skin tissue using PACO47174 at dilution of 1:100.