

CYREN Antibody, Biotin conjugated



PACO47177

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May act as a regulator of proteasome.

Gene ID:

CYREN

Uniprot

Q9BWK5

Synonyms:

Cell cycle regulator of non-homologous end joining (Cell cycle regulator of NHEJ)
(Modulator of retrovirus infection homolog), CYREN, C7orf49 MRI

Immunogen:

Recombinant Human Cell cycle regulator of non-homologous end joining protein (1-157AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A