

# ACTR5 Antibody, HRP conjugated



PACO47191

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Proposed core component of the chromatin remodeling INO80 complex which is involved in transcriptional regulation, DNA replication and probably DNA repair. Involved in DNA double-strand break repair and UV-damage excision repair.

**Gene ID:**

ACTR5

**Uniprot**

Q9H9F9

**Synonyms:**

Actin-related protein 5 (hARP5) (Sarcoma antigen NY-SAR-16), ACTR5, ARP5

**Immunogen:**

Recombinant Human Actin-related protein 5 protein (183-285AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A