

# OVOL2 Antibody, FITC conjugated



PACO47208

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Zinc-finger transcription repressor factor. Plays a critical role to maintain the identity of epithelial lineages by suppressing epithelial-to mesenchymal transition mainly through the up-regulation of ZEB1 expression. Positively regulates neuronal differentiation. Suppresses cell cycling and terminal differentiation of keratinocytes by directly repressing MYC and NOTCH1.

**Gene ID:**

OVOL2

**Uniprot**

Q9BRP0

**Synonyms:**

Transcription factor Ovo-like 2 (hOvo2) (Zinc finger protein 339), OVOL2, ZNF339

**Immunogen:**

Recombinant Human Transcription factor Ovo-like 2 protein (19-114AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A