

AKR1B1 Antibody, FITC conjugated



PACO47296

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the NADPH-dependent reduction of a wide variety of carbonyl-containing compounds to their corresponding alcohols with a broad range of catalytic efficiencies.

Gene ID:

AKR1B1

Uniprot

P15121

Synonyms:

Aldose reductase (AR) (EC 1.1.1.21) (Aldehyde reductase) (Aldo-keto reductase family 1 member B1), AKR1B1, ALDR1

Immunogen:

Recombinant Human Aldose reductase protein (114-230AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A