

AMACR Antibody, Biotin conjugated



PACO47301

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Racemization of 2-methyl-branched fatty acid, CoA esters. Responsible for the conversion of pristanoyl-CoA and C27-bile acyl-CoAs to their (S)-stereoisomers.

Gene ID:

AMACR

Uniprot

Q9UHK6

Synonyms:

Alpha-methylacyl-CoA racemase (EC 5.1.99.4) (2-methylacyl-CoA racemase), AMACR

Immunogen:

Recombinant Human Alpha-methylacyl-CoA racemase protein (13-175AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A