

ASS1 Antibody, Biotin conjugated



PACO47313

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Is indirectly involved in the control of blood pressure.

Gene ID:

ASS1

Uniprot

P00966

Synonyms:

Argininosuccinate synthase (EC 6.3.4.5) (Citrulline--aspartate ligase), ASS1, ASS

Immunogen:

Recombinant Human Argininosuccinate synthase protein (168-222AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A