

# DDX6 Antibody, FITC conjugated



PACO47376

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

In the process of mRNA degradation, plays a role in mRNA decapping. Blocks autophagy in nutrient-rich conditions by repressing the expression of ATG-related genes through degradation of their transcripts.

**Gene ID:**

DDX6

**Uniprot**

P26196

**Synonyms:**

Probable ATP-dependent RNA helicase DDX6 (EC 3.6.4.13) (ATP-dependent RNA helicase p54) (DEAD box protein 6) (Oncogene RCK), DDX6, HLR2 RCK

**Immunogen:**

Recombinant Human Probable ATP-dependent RNA helicase DDX6 protein (1-149AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A