

EPHB6 Antibody, Biotin conjugated



PACO47405

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Kinase-defective receptor for members of the ephrin-B family. Binds to ephrin-B1 and ephrin-B2. Modulates cell adhesion and migration by exerting both positive and negative effects upon stimulation with ephrin-B2. Inhibits JNK activation, T-cell receptor-induced IL-2 secretion and CD25 expression upon stimulation with ephrin-B2.

Gene ID:

EPHB6

Uniprot

O15197

Synonyms:

Ephrin type-B receptor 6 (HEP) (Tyrosine-protein kinase-defective receptor EPH-6), EPHB6

Immunogen:

Recombinant Human Ephrin type-B receptor 6 protein (410-553AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A