

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Plays an important role in the regulation of embryonic development, cell proliferation, cell differentiation and cell migration. Required for normal brain, eye, ear and limb development during embryogenesis. Required for normal development of the gonadotropin-releasing hormone (GnRH) neuronal system. Plays a role in neurite outgrowth in hippocampal cells.

Gene ID:

FGF8

Uniprot

P55075

Synonyms:

Fibroblast growth factor 8 (FGF-8) (Androgen-induced growth factor) (AIGF) (Heparin-binding growth factor 8) (HBGF-8), FGF8, AIGF


Immunogen:

Recombinant Human Fibroblast growth factor 8 protein (23-143AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human ovarian cancer using PACO47410 at dilution of 1:100.