

GSTM3 Antibody, HRP conjugated



PACO47443

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles. May govern uptake and detoxification of both endogenous compounds and xenobiotics at the testis and brain blood barriers.

Gene ID:

GSTM3

Uniprot

P21266

Synonyms:

Glutathione S-transferase Mu 3 (EC 2.5.1.18) (GST class-mu 3) (GSTM3-3) (hGSTM3-3), GSTM3, GST5

Immunogen:

Recombinant Human Glutathione S-transferase Mu 3 protein (1-225AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A