

# GUCY2C Antibody, HRP conjugated



PACO47447

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Receptor for the E.coli heat-stable enterotoxin (E.coli enterotoxin markedly stimulates the accumulation of cGMP in mammalian cells expressing GC-C). Also activated by the endogenous peptides guanylin and uroguanylin.

**Gene ID:**

GUCY2C

**Uniprot**

P25092

**Synonyms:**

Heat-stable enterotoxin receptor (STA receptor) (hSTAR) (EC 4.6.1.2) (Guanylyl cyclase C) (GC-C) (Intestinal guanylate cyclase), GUCY2C, GUC2C STAR

**Immunogen:**

Recombinant Human Heat-stable enterotoxin receptor protein (636-891AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A