

HRAS Antibody, Biotin conjugated



PACO47469

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Ras proteins bind GDP/GTP and possess intrinsic GTPase activity.

Gene ID:

HRAS

Uniprot

P01112

Synonyms:

GTPase HRas (H-Ras-1) (Ha-Ras) (Transforming protein p21) (c-H-ras) (p21ras) [Cleaved into: GTPase HRas, N-terminally processed], HRAS, HRAS1

Immunogen:

Recombinant Human GTPase HRas protein (12-181AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A