

PITX3 Antibody, HRP conjugated



PACO47475

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

PITX3

Uniprot

O75364

Synonyms:

Pituitary homeobox 3 (Homeobox protein PITX3) (Paired-like homeodomain transcription factor 3), PITX3, PTX3

Immunogen:

Recombinant Human Pituitary homeobox 3 protein (10-200AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A