

LPAR1 Antibody, Biotin conjugated



PACO47513

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

LPAR1

Uniprot

Q92633

Synonyms:

Lysophosphatidic acid, receptor 1 (LPA receptor 1) (LPA-1) (Lysophosphatidic acid, receptor Edg-2), LPAR1, EDG2 LPA1

Immunogen:

Recombinant Human Lysophosphatidic acid, receptor 1 protein (316-364AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A