

MAPK10 Antibody, Biotin conjugated



PACO47541

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

MAPK10

Uniprot

P53779

Synonyms:

Mitogen-activated protein kinase 10 (MAP kinase 10) (MAPK 10) (EC 2.7.11.24) (MAP kinase p49 3F12) (Stress-activated protein kinase 1b) (SAPK1b) (Stress-activated protein kinase JNK3) (c-Jun N-terminal kinase 3), MAPK10, JNK3 JNK3A PRKM10 SAPK1B

Immunogen:

Recombinant Human Mitogen-activated protein kinase 10 protein (133-316AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A