

# MTTP Antibody, Biotin conjugated



PACO47565

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Catalyzes the transport of triglyceride, cholesteryl ester, and phospholipid between phospholipid surfaces. Required for the secretion of plasma lipoproteins that contain apolipoprotein B.

**Gene ID:**

MTTP

**Uniprot**

P55157

**Synonyms:**

Microsomal triglyceride transfer protein large subunit, MTTP, MTP

**Immunogen:**

Recombinant Human Microsomal triglyceride transfer protein large subunit protein (19-151AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A