

PLA2G4A Antibody, HRP conjugated



PACO47599

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Selectively hydrolyzes arachidonyl phospholipids in the sn-2 position releasing arachidonic acid, Together with its lysophospholipid activity, it is implicated in the initiation of the inflammatory response.

Gene ID:

PLA2G4A

Uniprot

P47712

Synonyms:

Cytosolic phospholipase A2 (cPLA2) (Phospholipase A2 group IVA) [Includes: Phospholipase A2 (EC 3.1.1.4) (Phosphatidylcholine 2-acylhydrolase); Lysophospholipase (EC 3.1.1.5)], PLA2G4A, CPLA2 PLA2G4

Immunogen:

Recombinant Human Cytosolic phospholipase A2 protein (379-533AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A