

RAB7A Antibody, HRP conjugated



PACO47615

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

RAB7A

Uniprot

P51149

Synonyms:

Ras-related protein Rab-7a, RAB7A, RAB7

Immunogen:

Recombinant Human Ras-related protein Rab-7a protein (110-177AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A