

RAB7A Antibody, Biotin conjugated



PACO47617

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	RAB7A
Source:	Uniprot
Rabbit	P51149
Isotype:	Synonyms:
IgG	Ras-related protein Rab-7a, RAB7A, RAB7
Applications:	Immunogen:
ELISA	Recombinant Human Ras-related protein Rab-7a protein (110-177AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A