

RPS6 Antibody, FITC conjugated



PACO47624

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May play an important role in controlling cell growth and proliferation through the selective translation of particular classes of mRNA.

Gene ID:

RPS6

Uniprot

P62753

Synonyms:

40S ribosomal protein S6 (Phosphoprotein NP33) (Small ribosomal subunit protein eS6), RPS6

Immunogen:

Recombinant Human 40S ribosomal protein S6 protein (35-229AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A