

SH3BP5 Antibody, FITC conjugated



PACO47652

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Inhibits the auto- and transphosphorylation activity of BTK. Plays a negative regulatory role in BTK-related cytoplasmic signaling in B-cells. May be involved in BCR-induced apoptotic cell death.

Gene ID:

SH3BP5

Uniprot

O60239

Synonyms:

SH3 domain-binding protein 5 (SH3BP-5) (SH3 domain-binding protein that preferentially associates with BTK), SH3BP5, SAB

Immunogen:

Recombinant Human SH3 domain-binding protein 5 protein (49-455AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A