

SNAI2 Antibody, Biotin conjugated



PACO47685

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

SNAI2

Uniprot

O43623

Synonyms:

Zinc finger protein SNAI2 (Neural crest transcription factor Slug) (Protein snail homolog 2), SNAI2, SLUG SLUGH

Immunogen:

Recombinant Human Zinc finger protein SNAI2 protein (12-124AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A