

TOR1A Antibody



PACO47726

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:**Gene ID:**

TOR1A

Uniprot

O14656

Synonyms:

Torsin-1A (Dystonia 1 protein) (Torsin ATPase-1A) (EC 3.6.4) (Torsin family 1 member A),
TOR1A, DQ2 DYT1 TA TORA


Immunogen:

Recombinant Human Torsin-1A protein (88-125AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO47726 at dilution of 1:100.



Immunofluorescent analysis of HepG2 cells using PACO47726 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).