

TOR1A Antibody, FITC conjugated



PACO47728

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

TOR1A

Uniprot

O14656

Synonyms:

Torsin-1A (Dystonia 1 protein) (Torsin ATPase-1A) (EC 3.6.4) (Torsin family 1 member A), TOR1A, DQ2 DYT1 TA TORA

Immunogen:

Recombinant Human Torsin-1A protein (88-125AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A