

UNG Antibody, HRP conjugated



PACO47739

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Excises uracil residues from the DNA which can arise as a result of misincorporation of dUMP residues by DNA polymerase or due to deamination of cytosine.

Gene ID:

UNG

Uniprot

P13051

Synonyms:

Uracil-DNA glycosylase (UDG) (EC 3.2.2.27), UNG, DGU UNG1 UNG15

Immunogen:

Recombinant Human Uracil-DNA glycosylase protein (114-224AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A