

ALDOC Antibody, HRP conjugated



PACO47771

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytoskeleton, cytosol, extracellular exosome, cytoskeletal protein binding, fructose-bisphosphate aldolase activity, canonical glycolysis, epithelial cell differentiation, fructose 1,6-bisphosphate metabolic process, fructose metabolic process, gluconeogenesis.

Gene ID:

ALDOC

Uniprot

P09972

Synonyms:

Fructose-bisphosphate aldolase C (EC 4.1.2.13) (Brain-type aldolase), ALDOC, ALDC

Immunogen:

Recombinant Human Fructose-bisphosphate aldolase C protein (190-354AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A