

PET100 Antibody, FITC conjugated



PACO47852

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

PET100

Uniprot

P0DJ07

Synonyms:

Protein PET100 homolog, mitochondrial, PET100, C19orf79

Immunogen:

Recombinant Human Protein PET100 homolog, mitochondrial protein (26-73AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A