

MEI4 Antibody, HRP conjugated



PACO47863

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Required for meiotic DNA double-strand break formation.

Gene ID:

MEI4

Uniprot

A8MW99

Synonyms:

Meiosis-specific protein MEI4, MEI4

Immunogen:

Recombinant Human Meiosis-specific protein MEI4-like protein (22-137AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A