

C1QTNF9B Antibody, HRP conjugated



PACO47867

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

extracellular exosome.

Gene ID:

C1QTNF9B

Uniprot

B2RNN3

Synonyms:

Complement C1q and tumor necrosis factor-related protein 9B (C1q/TNF-related protein 9B) (CTRP9B) (Complement C1q and tumor necrosis factor-related protein 9-like), C1QTNF9B

Immunogen:

Recombinant Human Complement C1q and tumor necrosis factor-related protein 9B protein (2-192AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A