

# KLF5 Antibody, FITC conjugated



PACO47876

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Transcription factor that binds to GC box promoter elements. Activates the transcription of these genes.

**Gene ID:**

KLF5

**Uniprot**

Q13887

**Synonyms:**

Krueppel-like factor 5 (Basic transcription element-binding protein 2) (BTE-binding protein 2) (Colon krueppel-like factor) (GC-box-binding protein 2) (Intestinal-enriched krueppel-like factor) (Transcription factor BTEB2), KLF5, BTEB2 CKLF IKLF

**Immunogen:**

Recombinant Human Krueppel-like factor 5 protein (69-346AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A