

SWI5 Antibody, HRP conjugated



PACO47915

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Component of the SWI5-SFR1 complex, a complex required for double-strand break repair via homologous recombination.

Gene ID:

SWI5

Uniprot

Q1ZZU3

Synonyms:

DNA repair protein SWI5 homolog (HBV DNAPTP1-transactivated protein A) (Protein SAE3 homolog), SWI5, C9orf119 SAE3

Immunogen:

Recombinant Human DNA repair protein SWI5 homolog protein (54-112AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A