

C2CD3 Antibody, Biotin conjugated



PACO47929

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Component of the centrioles that acts as a positive regulator of centriole elongation. Promotes assembly of centriolar distal appendage, a structure at the distal end of the mother centriole that acts as an anchor of the cilium, and is required for recruitment of centriolar distal appendages proteins CEP83, SCLT1, CEP89, FBF1 and CEP164. Not required for centriolar satellite integrity or RAB8 activation. Required for primary cilium formation. Required for sonic hedgehog/SHH signaling and for proteolytic processing of GLI3.

Gene ID:

C2CD3

Uniprot

Q4AC94

Synonyms:

C2 domain-containing protein 3, C2CD3

Immunogen:

Recombinant Human C2 domain-containing protein 3 protein (2079-2265AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A