

TM4SF20 Antibody



PACO47938

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

integral component of membrane.

Gene ID:

TM4SF20

Uniprot

Q53R12

Synonyms:

Transmembrane 4 L6 family member 20, TM4SF20


Immunogen:

Recombinant Human Transmembrane 4 L6 family member 20 protein (114-184AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO47938 at dilution of 1:100.



Immunofluorescent analysis of HepG2 cells using PACO47938 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).