

TM4SF20 Antibody, Biotin conjugated



PACO47941

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

integral component of membrane.

Gene ID:

TM4SF20

Uniprot

Q53R12

Synonyms:

Transmembrane 4 L6 family member 20, TM4SF20

Immunogen:

Recombinant Human Transmembrane 4 L6 family member 20 protein (114-184AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A