

PACO47958

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

defense response.

**Gene ID:**

HCP5

**Uniprot**

Q6MZN7

**Synonyms:**

HLA class I histocompatibility antigen protein P5 (HLA complex protein P5) (Protein P5-1), HCP5

**Immunogen:**


Recombinant Human HLA class I histocompatibility antigen protein P5 protein (8-101AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human spleen tissue using PACO47958 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human lymph node tissue using PACO47958 at dilution of 1:100.