

PKD1L3 Antibody, Biotin conjugated



PACO47977

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Component of a calcium channel. May act as a sour taste receptor by forming a calcium channel with PKD1L3 in gustatory cells; however, its contribution to sour taste perception is unclear in vivo and may be indirect.

Gene ID:

PKD1L3

Uniprot

Q7Z443

Synonyms:

Polycystic kidney disease protein 1-like 3 (PC1-like 3 protein) (Polycystin-1L3), PKD1L3

Immunogen:

Recombinant Human Polycystic kidney disease protein 1-like 3 protein (22-194AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A