

SLC51B Antibody, HRP conjugated



PACO47983

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Essential component of the Ost-alpha/Ost-beta complex, a heterodimer that acts as the intestinal basolateral transporter responsible for bile acid, export from enterocytes into portal blood. Efficiently transports the major species of bile acid, . Modulates SLC51A glycosylation, membrane trafficking and stability activities.

Gene ID:

SLC51B

Uniprot

Q86UW2

Synonyms:

Organic solute transporter subunit beta (OST-beta) (Solute carrier family 51 subunit beta), SLC51B, OSTB

Immunogen:

Recombinant Human Organic solute transporter subunit β protein (57-128AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A