

LRRTM4 Antibody, FITC conjugated



PACO47988

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in the development and maintenance of the vertebrate nervous system. Exhibits strong synaptogenic activity, restricted to excitatory presynaptic differentiation.

Gene ID:

LRRTM4

Uniprot

Q86VH4

Synonyms:

Leucine-rich repeat transmembrane neuronal protein 4, LRRTM4

Immunogen:

Recombinant Human Leucine-rich repeat transmembrane neuronal protein 4 protein (259-419AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A