

COMMD1 Antibody, HRP conjugated



PACO48007

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

COMMD1

Uniprot

Q8N668

Synonyms:

COMM domain-containing protein 1 (Protein Murr1), COMMD1, C2orf5 MURR1

Immunogen:

Recombinant Human COMM domain-containing protein 1 protein (2-190AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A