

ZDHHC1 Antibody, HRP conjugated



PACO48023

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

endoplasmic reticulum, extracellular exosome, DNA binding, palmitoyltransferase activity, protein palmitoylation.

Gene ID:

ZDHHC1

Uniprot

Q8WTX9

Synonyms:

Probable palmitoyltransferase ZDHHC1 (EC 2.3.1.225) (DHHC domain-containing cysteine-rich protein 1) (Zinc finger DHHC domain-containing protein 1) (DHHC-1) (Zinc finger protein 377), ZDHHC1, C16orf1 ZNF377

Immunogen:

Recombinant Human Probable palmitoyltransferase ZDHHC1 protein (283-389AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A