

RPTOR Antibody



PACO48030

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Gene ID:

RPTOR

Uniprot

Q8N122

Synonyms:

Regulatory-associated protein of mTOR (Raptor) (p150 target of rapamycin (TOR)-scaffold protein), RPTOR, KIAA1303 RAPTOR


Immunogen:

Recombinant Human Regulatory-associated protein of mTOR protein (13-128AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO48030 at dilution of 1:100.



Immunofluorescent analysis of Hela cells using PACO48030 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO48030 at dilution of 1:100.