

UNC13C Antibody, HRP conjugated



PACO48035

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May play a role in vesicle maturation during exocytosis as a target of the diacylglycerol second messenger pathway. May be involved in the regulation of synaptic transmission at parallel fiber - Purkinje cell synapses.

Gene ID:

UNC13C

Uniprot

Q8NB66

Synonyms:

Protein unc-13 homolog C (Munc13-3), UNC13C

Immunogen:

Recombinant Human Protein unc-13 homolog C protein (414-645AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A