

CEP192 Antibody, FITC conjugated



PACO48052

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Required for mitotic centrosome and spindle assembly. Appears to be a major regulator of pericentriolar material (PCM) recruitment, centrosome maturation, and centriole duplication.

Gene ID:

CEP192

Uniprot

Q8TEP8

Synonyms:

Centrosomal protein of 192 kDa (Cep192), CEP192, KIAA1569

Immunogen:

Recombinant Human Centrosomal protein of 192 kDa protein (447-711AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A