

DDIT4L Antibody, Biotin conjugated



PACO48061

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Inhibits cell growth by regulating the TOR signaling pathway upstream of the TSC1-TSC2 complex and downstream of AKT1.

Gene ID:

DDIT4L

Uniprot

Q96D03

Synonyms:

DNA damage-inducible transcript 4-like protein (HIF-1 responsive protein RTP801-like) (Protein regulated in development and DNA damage response 2) (REDD-2), DDIT4L, REDD2 RTP801L

Immunogen:

Recombinant Human DNA damage-inducible transcript 4-like protein (1-163AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A