

ATXN2 Antibody, Biotin conjugated



PACO48085

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in EGFR trafficking, acting as negative regulator of endocytic EGFR internalization at the plasma membrane.

Gene ID:

ATXN2

Uniprot

Q99700

Synonyms:

Ataxin-2 (Spinocerebellar ataxia type 2 protein) (Trinucleotide repeat-containing gene 13 protein), ATXN2, ATX2 SCA2 TNRC13

Immunogen:

Recombinant Human Ataxin-2 protein (174-383AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A