

PPP1R12C Antibody, Biotin conjugated



PACO48117

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Regulates myosin phosphatase activity.

Gene ID:

PPP1R12C

Uniprot

Q9BZL4

Synonyms:

Protein phosphatase 1 regulatory subunit 12C (Protein phosphatase 1 myosin-binding subunit of 85 kDa) (Protein phosphatase 1 myosin-binding subunit p85), PPP1R12C, LENG3 MBS85

Immunogen:

Recombinant Human Protein phosphatase 1 regulatory subunit 12C protein (317-450AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A