

PIM2 Antibody, HRP conjugated



PACO48127

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

PIM2

Uniprot

Q9P1W9

Synonyms:

Serine/threonine-protein kinase pim-2 (EC 2.7.11.1) (Pim-2h), PIM2

Immunogen:

Recombinant Human Serine/threonine-protein kinase pim-2 protein (1-210AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A