

LZTFL1 Antibody, HRP conjugated



PACO48159

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Regulates ciliary localization of the BBSome complex. Together with the BBSome complex, controls SMO ciliary trafficking and contributes to the sonic hedgehog (SHH) pathway regulation. May play a role in neurite outgrowth. May have tumor suppressor function.

Gene ID:

LZTFL1

Uniprot

Q9NQ48

Synonyms:

Leucine zipper transcription factor-like protein 1, LZTFL1

Immunogen:

Recombinant Human Leucine zipper transcription factor-like protein 1 protein (116-299AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A