

RAD54B Antibody, HRP conjugated



PACO48183

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in DNA repair and mitotic recombination. May play an active role in recombination processes in concert with other members of the RAD52 epistasis group.

Gene ID:

RAD54B

Uniprot

Q9Y620

Synonyms:

DNA repair and recombination protein RAD54B (EC 3.6.4) (RAD54 homolog B), RAD54B

Immunogen:

Recombinant Human DNA repair and recombination protein RAD54B protein (1-167AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A